

Study Cases of PackGene's AAV Services

Case 1

Products: AAV8 and AAV9 from PackGene Biotech

Injection Method: 300nl was injected bilaterally into the anterior cingulate cortex

Recommended Dose: 50ul, 5×10^{11} GC/ml


Targeting Site: Anterior Cingulate Cortex

Animal Model: Global Shank3 KO and Conditional ACC Shank3 KO Mouse Model

Journal: Nature Neuroscience, 2022 (IF=28.771)

Paper Title: Anterior cingulate cortex dysfunction underlies social deficits in Shank3 mutant mice

DOI: <https://doi.org/10.1038/s41593-019-0445-9>



SaCas9 plasmid construction and AAV vector production.

To minimize offtarget effects, three guide RNAs (gRNA1: GGGCTATTCCAGCCTCCCTCC; gRNA2: GGCAGGGGCGTGTCCAGGTTAG; gRNA3: TTGGCGGCCACACG GGCGCGG) corresponding to SHANK3 were designed using the CRISPR design tool (<http://www.rgenome.net/cas-designer/>) and cloned into pX601-AAVCMV::NLS-SaCas9-NLS-3xHA-bGHpA;U6::Bsal-sgRNA... After the plasmids were confirmed by sequencing, we transfected those constructs into N2A mouse cells. A T7E1 assay was performed to evaluate indel efficiency at 72h after transduction. We then chose gRNA2, which had the best efficiency, to package in AAV serotype 8 (PackGene Biotech) for administration to the mice.

Case 2

Products: HBV8-B-AAV from PackGene Biotech

Injection Method: Tail Vein Injection


Targeting Site: Liver

Animal Model: C57BL/6

Journal: Nature Neuroscience, 2022 (IF=10.103)

Paper Title: A broad-spectrum nanobody targeting the C-terminus of the hepatitis B surface antigen for chronic hepatitis B infection therapy

DOI: <https://doi.org/10.1016/j.antiviral.2022.105265>



All mouse experiments were performed according to the guidelines of the Laboratory Animal Management Ethics Committee at Xiamen University. C57BL/6 mice were injected with HBV/adeno-associated virus (AAV) (HBV-B from PackGene Biotech, HBV-D from Five Plus Molecular Medicine Institute, China) via the caudal vein. HCAb treatment was conducted as previously described (Zhou et al., 2020).

Case 3

Products: AAV-PhP.eB from PackGene Biotech

Injection Method: Retro-ocular Injection

Recommended Dose: 1.3×10^{12} GC/ml, 100 μ l


Targeting Site: Neuron

Animal Model: Shank3B KO mice; Sst-IRES-Cre knock-in mice

Journal: Frontiers in cellular neuroscience, 2020 (IF=6.147)

Paper Title: A whole-brain cell-type-specific sparse neuron labeling method and its application in a shank3 autistic mouse model

DOI: <https://doi.org/10.3389/fncel.2020.00145>



Retro-Ocular Injection

Mice were given retro-ocular injection of AAV as previously described. Briefly, mice were anesthetized with 3% isoflurane. Then, 100 μ l of AAV php.eb-DIO-EGFP-P2A-EGFPf with a titer of 1.3×10^{12} gc/ml (PackGene Biotech, LLC, China) was injected into the retro-orbital sinus with a 27-G needle and a 1-ml syringe. The mice were then placed in

a warm and moist environment to wait for resuscitation. After resuscitation, the mice were put back in the cage and housed for 3 weeks before the next step.

Case 4

Products: AAV8 from PackGene Biotech

Injection Method: Ear Vein Injection

Recommended Dose: AAV8-SpCas9 3×10^{12} GC; AAV8-sgRNA4-Donor 6×10^{12} GC


Targeting Site: Liver

Animal Model: Rabbits

Journal: Molecular Therapy, 2021 (IF=12.910)

Paper Title: CRISPR/Cas9-Mediated Gene Correction in New born Rabbits with Hereditary Tyrosinemia Type I

DOI: <https://doi.org/10.1016/j.ymthe.2020.11.023>



AAV8 Production and Purification

AAV8 was produced and purified by PackGene Biotech. Donor template and optimal sgRNA with U6 promoter were packaged into an AAV8 with mCherry fluorescent reporter, referred to as AAV8- sgRNA4-Donor. SpCas9 with the miniCMV promoter was packaged into the other AAV8, referred to as AAV8- SpCas9. The purified AAV8 was stored in a freezer at -80°C .

Case 5

Products: AAVrh10 from PackGene Biotech

Injection Method: Intranasal Administration


Recommended Dose: 2×10^{11} GC/ml, 40 μ l

Animal Model: 6-8-week-old female BALB/c mice

Journal: Antiviral Research, 2022 (IF=10.103)

Paper Title: Rational design of AAVrh10-vectored ACE2 functional domain to broadly block the cell entry of SARS-CoV-2 variants

DOI: <https://doi.org/10.1016/j.antiviral.2022.105383>



AAVrh10 expressing the different ACE2 fragments, including shACE2, shACE2PD, shACE2(PD + Neck), hACE224-83, dACE224-83, were packaged by PackGene Biotech.

Case 6

Products: AAV-PHP.eB from PackGene Biotech

Injection Method: Cochlear Injection

Recommended Dose: 1.43×10^{10} GC (2.86×10^{13} GC/ml in 500nl/ 3.99×10^{13} GC/ml in 358.4nl)


Targeting Site: Cochlea

Animal Model: Myo6^{C442Y/+} and Myo6^{C442Y/C442Y} mice

Journal: Antiviral Research, 2022 (IF=10.103)

Paper Title: Rescue of autosomal dominant hearing loss by in vivo delivery of mini dCas13X-derived RNA base editor

DOI: 10.1126/scitranslmed.abn0449



AAV virus production

AAV-PHP.eB serotype (25, 35) was used in this study. The mxABE plasmid with target or nontarget gRNA was sequenced before packaging into AAV-PHP.eB vehicle, and the AAV vectors were packaged by PackGene Biotech. The vector titer was 2.86×10^{13} and 3.99×10^{13} genome copies/ml as determined by qPCR specific for the inverted terminal repeat of the AAV-PHP.eB-Myo6 and control virus.

Case 7

Products: AAV9 from PackGene Biotech

Injection Method: Intratracheal Injection

Recommended Dose: 1×10^{10} vg/ml


Targeting Site: Lung

Animal Model: C57BL/6 BMP4^{+/+} and BMP4^{+/-} mice

Journal: European respiratory journal, 2022 (IF=33.795)

Paper Title: Bone morphogenetic protein 4 inhibits pulmonary fibrosis by modulating cellular senescence and mitophagy in lung fibroblasts

DOI: 10.1183/13993003.02307-2021



heterozygous null (BMP4^{+/-}) mice for all the experiments. In addition, separate batches of mice were used for the prevention and therapeutic studies. In the prevention study, mice were administrated with BMP4-expressing AAV9 viral genome particle (1.0×10^{10} vg/mL, PackGene Biotech, Guangzhou, China) via intratracheal injections. After 21 days post-viral treatment, these mice were...

Case 8

Products: CRISPR AAV vectors from PackGene Biotech

Injection Method: Intramuscular Injection

Recommended Dose: 25 μ L, 1.5E+12 vg/ml


Targeting Site: muscle

Animal Model: DMD mouse model

Journal: Genome Medicine, 2021 (IF=15.266)

Paper Title: In vivo genome editing in mouse restores dystrophin expression in Duchenne muscular dystrophy patient muscle fibers

DOI: 10.1186/s13073-021-00876-0



AAV production and injection

CRISPR AAV vectors were generated by PackGene Biotech Co. (Guangzhou, China). For intramuscular injection into NSI mice, the animals were anesthetized by abdominal injection of 5% chloral hydrate, followed by injection of CRISPR AAV vector (25 μ L, 1.5E+12 vg) into the TA muscle transplanted with DMD–MDSCs.

Case 9

Products: AAV9 and AAV5 from PackGene Biotech

Injection Method: Intramuscular Injection

Recommended Dose: 2×10^{13} GC/ml, 1 μ l


Targeting Site: Central Nervous System

Animal Model: Adult C57BL/6J Mice

Journal: Neural Regeneration Research, 2021 (IF=6.058)

Paper Title: Lineage tracing of direct astrocyte-to-neuron conversion in the mouse cortex

DOI: 10.4103/1673-5374.295925



AAV serotype 9 (AAV9) and 5 (AAV5) were produced by PackGene Biotech, LLC, purified through iodixanol gradient ultracentrifuge and subsequent concentration.

Case 10


Products: AAV2 from PackGene Biotech

Animal Model: Hepatitis B Model Mice

Journal: Frontiers in Immunology, 2022 (IF=8.786)

Paper Title: PreS/2-21-Guided siRNA Nanoparticles Target to Inhibit Hepatitis B Virus Infection and Replication

DOI: 10.3389/fimmu.2022.856463



The recombinant adeno-associated virus AAV-HBV-002 applied to construct the hepatitis B model was purchased from PackGene Biotech, Co., Ltd. AAV-HBV-002 containing 1.3× HBV genome, is characterized by the production of HBV DNA, HBeAg and HBsAg, and the genotype is C2 and serotype is adr. Sixteen mice

Case 11

Products: AAV9 from PackGene Biotech

Injection Method: Tail Vein Injection

Recommended Dose: 1×10^{11} vg/mouse


Targeting Site: Heart

Animal Model: Adult Male Kunming Mice (weight 22–25g)

Journal: Molecular Therapy Nucleic Acids, 2020 (IF=10.183)

Paper Title: MicroRNA-20b Promotes Cardiac Hypertrophy by the Inhibition of Mitofusin 2-Mediated Inter-organelle Ca^{2+} Cross-Talk

DOI: <https://doi.org/10.1016/j.omtn.2020.01.017>



Synthesis and Administration of rAAV9-anti-miR-20bAAV was used as carrier, among which rAAV9 is the most efficient vector for myocardial transduction. For this regard, rAAV9-antimiR-20b or rAAV9-NC was produced (PackGene Biotech, Guangzhou, China) and the sequences were delivered into mouse, respectively, through tail vein injection at 1×10^{11} vg (viral genomes) per animal at 5 weeks before TAC. Animals were divided into four

Case 12

Products: All AAV from PackGene Biotech


Targeting Site: Pancreas

Animal Model: CP Model Mice

Journal: Journal of Translational Medicine, 2022 (IF=8.440)

Paper Title: St13 protects against disordered acinar cell arachidonic acid pathway in chronic pancreatitis

DOI: 10.1186/s12967-022-03413-8



Blank adenovirus was used as negative control. All adeno-associated virus (PackGene Biotech, China) used are listed in Additional file 4: Fig. S4. Please see the Supplementary materials and methods for more details.